



# Meet Your New PERS

## Digital Emergency Response System

The Digital PERS device was created to provide you comfort, reassurance and peace of mind by making it possible for help to be just a touch of a button away.

### Get to Know Your Unit



- 1 Speaker**  
Used to listen to the Care Center
- 2 Emergency Button**  
Press to send emergency alarms to the Care Center
- 3 Alarm Status Light**  
Dim white light indicates normal operation  
Bright white light indicates calling the Care Center  
Flashing bright white light indicates unit is in programming mode
- 4 Test Button**  
Used to send test signal to the Care Center
- 5 Reset Button**  
Used to cancel alarm
- 6 Power Light**  
Solid light indicates unit is properly connected  
Flashing light indicates AC power is not available and running on its back-up battery
- 7 Trouble Light**  
Indicates phone line trouble when flashing
- 8 Microphone**  
Used to communicate with the Care Center
- 9 On/Off Button**  
Switch on to power the unit
- 10 Power Jack**  
Connect with the provided AC adapter to a power outlet
- 11 Phone Jack**  
Connects your unit to your telephone
- 12 Wall Jack**  
Connects the unit to your telephone wall jack or modem
- 13 Programming Button**  
Used to range test and learn in emergency sensors

### Unit Facts & Features:

- ◆ Fall detection available
- ◆ Water-resistant button
- ◆ The button can be worn around your wrist or neck
- ◆ Please test your unit monthly to ensure it is working properly
- ◆ Unit will perform silent self tests to the Care Center
- ◆ For additional questions, call tech support at [1-800-860-4230](tel:1-800-860-4230)