


Meet Your New PERS

MXD Emergency Response System

MXD was created to provide you comfort, reassurance & peace of mind by making it possible for help to be just a touch of a button away. The MXD was designed to be easy to install and easy to maintain.

Get to Know Your Unit



- 1 Microphone**
Used to talk to the Care Center
- 2 Help Button**
Press to call the Care Center (works just like the button)
- 3 Speaker**
Used to hear the Care Center
- 4 Status Light**
Displays the current operating condition
- 5 Wall Jack**
Used to connect the unit to your home phone service
- 6 Phone Jack**
For service provider use only
- 7 Reset Button**
Used to Reset the unit
- 8 T/L Button**
The Test and Learn Button. Used in the programming and testing of the unit

Unit Facts & Features:

- ✦ The button can be worn around your neck or as a wristband
- ✦ The button has a range of up to 600 feet and an average 3-5 year battery life
- ✦ Compatible with traditional phone lines, cable phones or VoIP systems.
- ✦ Care Center monitors unit's battery status
- ✦ Water-resistant button
- ✦ 60+ hours of backup battery in case of AC power loss
- ✦ Optional fall detection pendant can be added
- ✦ Optional Medication monitoring can be added
- ✦ Status Light Colors:
 - **Solid Green-** Unit is properly connected
 - **Flashing Green-** AC power is not available, unit is running on back up battery
- **Steady Red-** Unit is trying to call the Care Center
- **Flashing Red-** Unit isn't connected to telephone service
- **Flashing Orange-** Unit is operating on backup battery and is not connected to working telephone service
- ✦ For additional questions, call tech support at **1-800-860-4230**