


Quick Install

MedReady Medication Dispenser System

Welcome to your new medication dispenser system. If you have any questions or in need of assistance, do not hesitate to call us at [1-800-860-4230](tel:1-800-860-4230).


Connections

1. Flip MedReady upside down
2. Attach back up battery\replace battery lid
3. Attach supplied AC adapter to bottom of MedReady
4. Plug into household power outlet. Make sure that the outlet is **not** operated by a light switch or power strip
5. For “MedReady Modem” plug supplied phone cord into bottom of MedReady and wall jack



Setting Times/Alarms

1. Set clock to current time by holding the “Time Set” button while pressing Hour, Minute, and AM-PM buttons
2. To set alarm times; Press and Hold #1 alarm (1st alarm of the day) then use the Hour, Minute, and AM-PM buttons Repeat step for alarms 2-4 as needed
3. When active there will be an “ON” light above each active alarm time 1-4
4. Slide the white bar to the center near the “X” until it locks in place
5. Press alarm test, after alarm goes off, reset white bar back to the “X” location (This will silence alarm)



Loading the Dispenser

1. Fill the medication wheel one dose per compartment
2. Load in a counter clockwise direction
3. 4 doses a day will hold 7 days of meds
4. Slide the white bar back to the center near the “X”
5. Close the lid, lock it, and verify the alarm lights are “ON” for any dispenses you have programmed

✦ Twist lid counter clockwise to remove

MedReady Install Tips

- MedReady will automatically set the correct time after 20 minutes
- MedReady Connect should be pre-paired to your PERS device
- Always observe AM-PM on the clock & alarm times
- MedReady display will flash “:” when running on back up battery