


# Meet your new mPERS

MobileCare

- 1 Cellular Indicator**  
Indicates cell service level  
(Red, Yellow, Green)
- 2 Battery Indicator**  
Indicates the battery level  
(Red, Yellow, Green)
- 3 Call Button**  
Press to send emergency alarms to  
the Care Center
- 4 Speaker**  
Used to hear the Care Center
- 5 Microphone**  
Used to talk to the Care Center
- 6 Charging Cradle**  
Where you will place your Mobile Care  
to charge
- 7 Gold Charger Contacts**  
Contact with charging station, allows unit  
to charge



1400 Commerce Center Drive  
Franklin, OH 45005

V03042019

If you know someone who could  
use our assistance, call us at  
**1-800-860-4230** or log on to  
**www.vricares.com** for more  
information on our products  
and services.

## MobileCare Instructions

mobile medical alert systems




# MobileCare

## Instruction Guide




**1** Connect the power cord into the back of the charging cradle, then plug the MobileCare charging cradle into a power outlet. Make sure the outlet is NOT operated by a light switch or power strip.



**2** Place the MobileCare in the charging cradle and wait for a welcome message. The diamond shaped battery light will turn green when fully charged.



**3** While leaving the MobileCare in the charging cradle, **press** your call button to activate your unit.



**4** The Care Center will talk to you from your MobileCare mPERS unit. Let us know you are a new client testing your unit!



You're going places. Take us with you.

MobileCare was created to provide you comfort, reassurance and peace of mind by making it possible for help to be just a touch of a button away.


It is important to **test your device monthly** to ensure device is working properly. Before you test your unit, make sure it is **completely charged**.

Fall detection is not a perfect science. If you are able, **always** press your help button in the case of an emergency.

If you have any questions or are in need of assistance, do not hesitate to call us at **1-800-860-4230**.

### MobileCare Light Guide

Signal Light	●	Strong	Battery Life	●	Good
	●	Moderate		●	Moderate
	●	No Signal		●	Charge

 **Multiple Flashing Lights** - Device is rebooting and attempting to connect to the network. When this occurs, place device on charging cradle.