


Meet Your New Device

ForaCare BP Monitoring System

ForaCare BP was created to provide you comfort, reassurance and peace of mind by allowing you to check your blood pressure in the comfort of your own home.

Get to Know Your System



- 1 Cellular Signal Indicator**
Indicates a wireless transmission. Off in standby mode & blue during transmission.
- 2 USB Indicator**
(Used for programming only)
- 3 Power/Battery Indicator**
Indicates power status of the device. Dark when charged, flashing green when charging, & red when low.
- 4 S Button**
Used to enter and confirm the meter settings
- 5 M Button**
Used to enter meter memory
- 6 Average Button**
Displays BP averages. (Normally not used)
- 7 Cell Signal Button**
Press and hold to test cellular signal strength
- 8 On/Off Button**
Press to start taking a reading
- 9 Speaker**
Depending on your meter's setting, your meter may have audio prompts

System Facts & Features:

- ✦ The ForaCare device is intended for single patient use only and should not be shared
- ✦ Always make sure you take your reading on the same arm
- ✦ AT&T Cellular Service
- ✦ Avoid caffeine, alcohol & tobacco for 30 minutes before taking a reading
- ✦ Sit down for at least 10 minutes before taking a reading
- ✦ After taking a reading, the ForaCare device will automatically upload your reading via unit's cellular connection
- ✦ After taking a reading, the ForaCare device will automatically power off
- ✦ The ForaCare BP meter has a rechargeable battery, but should be plugged into a power outlet under normal use
- ✦ For additional questions, call tech support at **1-800-860-4230**