

## Glucometer and Polytel Device (if applicable)

- 1** Unpack the Polytel and Glucometer.



- 2** On the back of the Polytel Device open the battery compartment. Place batteries in the device and close the compartment.




- 3** Take your Glucometer reading.



- 4** Connect Polytel Device to the glucometer by plugging the cable into the port on the Glucometer.



- 5** When plugged in, press the button shown below to send your reading. Your glucometer will say PC.



- 6** Once you have set up all units applicable to you, please call our Vitals Monitoring Team at [1-888-875-0976](tel:1-888-875-0976) then press "1" to make sure your readings were received.